

Safety through precision

GLOOR
MEDICAL

SWISS optio terminal unit AGSS



terminal unit AGSS

The terminal unit AGSS at a glance:

- Attractive, functional design
- Full metal and robust
- Service-friendly and long-lasting
- For use in hospitals, clinics and emergency services

The **terminal unit AGSS** (**A**naesthetic **G**as **S**cavenging **S**ystem) serves to channel off in a controlled manner excess anaesthetic gas.

The negative pressure (vacuum), which is required by the **AGSS terminal unit** for the extraction of the excess and exhaled anaesthetic gas, is generated via an integrated ejector operated with medical air. The extraction is carried out using a hose which is connected via a connector (1L) corresponding to the EN ISO 9170-2 standard to the **AGSS terminal unit**.

By coupling the hose, the **AGSS terminal unit** will be put into operation. An optical operating indicator (green light) indicates operation of the **AGSS terminal unit**.


The **terminal unit AGSS** is for assembling in ducts or for flush mounted installations. Because it is made of solid metal, the **AGSS terminal unit** is very robust and correspondingly durable.

The **terminal unit AGSS** complies with the requirements of the standard EN ISO 9170-2. CE marking according to the Guideline for Medical Products 93/42/EEC.

Technical data :

Pressure range for the supply of the ejector	max. 6 bar
Model variants	in-duct assembly / ceiling supply units flush mounting
Operating gas ejector	air (medical)
Suction medium	anaesthesia gas
Connections	discharge suction medium : pipe end with outside diameter of 15 mm for connection with a flexible hose to the clinic's exhaust line supply of ejector : Ø4 mm pneumatic connector for the ejector's supply air
Material	brass (outer parts nickel plated), copper, POM, EPDM, hostaphan
Operating temperature	-20° to +60°C
Dimensions	W x D x diameter : 80 mm x 69.5 mm x Ø 48 mm (base assembly group)
Weight	900 gr

Model variants GM31xx



Necessary information when ordering :

Model variants in-duct assembly, flush-mounting

Available accessories :

Art. 661690 angle probe AGSS 1 litre Ø 22 mm

Subject to changes. All trademarks are property of Gloor Ltd.