


SWISS optio cylinder trolley



The **cylinder trolley** ensures the safe and convenient transport of single cylinders. The five solid wheels guarantee for reliable stability and easy maneuverability.

With integrated DIN rail for the customised attachment of equipment.

The high quality workmanship as well as the use of top quality components ensures that the **cylinder trolley** is very robust to match the high demands of everyday clinical usage.

Mainly intended for use in clinical departments and emergency units.

cylinder trolley

Technical data : GX10047/4

Material	basic frame aluminium, anodized natural, wheels Ø 75 mm, fixed castor INOX, bearing electrically conducting
Mounting rail	standard rail DIN25 / 10 mm (length 200 mm)
Cylinder holder	for 5 and 10 l cylinder (Ø 140 mm)
Dimensions	height : 970 mm, width : 515 mm, depth : 515 mm
Weight	4 kg (without cylinder)

Subject to changes. All trademarks are the property of Gloor Ltd.